

Application Guide

Shaped Pad style products (lying position)



1. Fold the pad lengthways with the backsheet turned outwards. Do not touch the inside of the pad.

2. Pass the folded pad from front to back, creating a channel shape between the legs. If necessary, place a towel between the knees to ease the process.



3. Pull out the back of the pad horizontally, and position the upper rear at the end of the coccyx.

4. Pull up the fitting pants.



5. With the wearer lying on their back fan out the pad on the front. The pad should create a channel shape in the groin area and fit snugly to the body.



6. Adjust the fitting pants at the front. If possible, the pant legs should be placed on the backsheet edges of the pad.

7. The wetness indicator is located on the reverse of the pad. And changes from yellow to blue when wet. Once the wetness indicator is around two thirds full, this can be used as a guide for when the product needs changing.

If there is a gap in the blue graduation scale on the wetness indicator then there could be a fitting problem.

Please note that in any cases of faecal incontinence, the pad should be changed as soon as possible.



Application Guide

Shaped Pad style products (standing position)



1. Fold the pad lengthways with the backsheet turned outwards. Do not touch the inside of the pad.
2. Ensuring the user is properly supported pass the pad from front to back. Position the upper rear end of the pad at the coccyx.

3. Fan out the back of the pad horizontally to create a channel shape between the legs.



4. Fan out pad horizontally at the front.



5. Pull the fitting pants into position. If possible, the pant legs should be placed on the backsheet edges of the pad.



6. Adjust the fitting pants at the front and make sure the edges of the pad ease into the groin area.

7. The wetness indicator is located on the reverse of the pad. And changes from yellow to blue when wet. Once the wetness indicator is around two thirds full, this can be used as a guide for when the product needs changing.

If there is a gap in the blue graduation scale on the wetness indicator then there could be a fitting problem.

Please note that in any cases of faecal incontinence, the pad should be changed as soon as possible.

