

TENA Men


Engineered especially for men

TENA Men, available in three absorbency levels, is especially designed to fit the male form – ensuring maximum security and discretion. The cup shaped pad follows the contours of the body to give a snug and secure fit. The adhesive strip secures the pad in position within normal underwear. The soft material is gentle to the skin and completely noise free. The unique Odour Control™ inhibits the development of unwanted smells for a more dignified way of life.

Anatomically shaped pad is designed to fit the male anatomy in order to provide maximum comfort and leakage security and reduce odours

Fast absorption, dry surface reduces skin irritation

Unique Odour Control means urine is locked away in the pad, preventing the development of unwanted odour for discretion and comfort




Single-wrapped products for added discretion and hygiene


Adhesive strip holds the pad in place in the underwear, minimising the risk of leakage

Leak proof breathable textile-like backing improves comfort and reduces skin irritation to prevent staining


TENA Men Fitting Guide



Open the TENA Men packaging and unfold the pad, removing the tab to expose the adhesive strip.



Make sure that the wider part of the pad is positioned upwards.



Place the pad in snug-fitting underwear so that the adhesive strip attaches to the inside front.



Do not use with loose-fitting underwear or boxer shorts.

Contact your TENA Account Manager or email hcmarketing@essity.com for more information about TENA Men.

	Level 1	Level 2	Level 3
Absorbency	●●●○	●●●●○	●●●●●
Product Code	750661	750797	750830
Working Absorbency	150ml	200ml	300ml

